TR

**ELMAS KAPLI DENTAL FREZ**  
KULLANMA KILAVUZU

ENG

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
## EN

### 1. PRODUCT INFORMATION

<b>Intended Use</b>	<p>The general purpose of these products is made of stainless steel coated with small diamond crystals. It is a rotary cutting device designed to fit into a dental handpiece and used to carve and cut hard structures in the mouth. It is used to cut hard tissue such as enamel and dentin.</p> <p>For products requiring detailed information, separate usage locations are available on our website <a href="http://www.bosphorus.com.tr">www.bosphorus.com.tr</a></p> <p>These details are provided on the general instructions.</p>
---------------------	--

<b>1</b>	<b>Candle Flame Dental Burs</b>	It is used for contour correction, shaving of the pointed tubercle (anatomical deformity of the tooth), composite leveling, cutting and correction.
<b>2</b>	<b>Cylinder Dental Burs</b>	It is used for cavity preparation, steps and corrections, abrasion of molars and cavity expansion
<b>3</b>	<b>Round Dental Bur</b>	<p>They are used to remove caries and when opening an endodontic or restoration cavity. It serves to make guide grooves on the tooth. It is also used in deep cavities to clean the colored but hard bruised area.</p> <p>If the cavity still continues after the cavity reaches a certain depth, then the cavity is not deepened any further, only the bruised part is cleaned locally.</p> <p>When a rond bur is used, a thinner and easily broken filling edge is formed.</p>
<b>4</b>	<b>Inverted Cone Dental Burs</b>	It is used for cavity preparation and cavity floor correction.
<b>5</b>	<b>Torpedo Dental Burs</b>	Torpedo burs are used according to the edge angle to be prepared in stepped tooth cutting (chamfer or shoulder). Knowing the diameter of the bur helps to determine the depth of the guide groove to be prepared and thus the thickness of the material to be reduced.
<b>6</b>	<b>Conical Dental Burs</b>	It is used in caries, cavity opening, cavity floor correction, tooth enamel and dentin cutting and correction.
<b>7</b>	<b>Needle Dental Burs</b>	Used for cutting tooth enamel and dentin, post-cutting corrections and opening cavities.
<b>8</b>	<b>Egg Dental Burs</b>	Composite correction, lingual side correction of front teeth. It is used to correct the labial back side and correct tooth enamel and dentin.
<b>9</b>	<b>Wheel Dental Burs</b>	It is used to create a guide groove, to prepare a guide groove on the cutting edge (labia in the front teeth) surface, and to prepare a guide groove on the front tooth surface.

**MILLING MODELS AND APPLICATIONS**

No	Product Model	Product Images	Length Measure	Colour											
					Prophylaxis	Crowns Preparations (Cutting)	Old Fill Removal (carving)	Crowns and Bridge Technique (Cutting)	Model Casting Technique (Cutting-Carving)	Cavity Preparations (Carving)	Mouth Surgery (Cutting)	Acrylic Technique (carving)	Tooth Orthopedics (Cutting)	Model Making (Cutting)	
1	Candle Flame		Standard Short L XL	Green Black Yellow Blue Red White	X	X		X	X						
2	Cylinder				X	X	X	X	X	X					
3	Round				X	X	X	X	X	X	X	X	X	X	
4	Inverted Cone				X		X	X	X	X					
5	Torpedo					X									
6	Conical				X	X	X				X				
7	Needle				X	X		X	X	X					X
8	Egg				X	X	X	X	X	X					
9	Wheel						X					X			

<b>Indications</b>	<p>Treatment of carious teeth</p> <p>Tooth preparation for prosthetic treatment</p> <p>Removal of fillings from already restored teeth</p> <p>Removal of denture</p> <p>Separating of teeth or denture</p>
<b>Contraindications</b>	<p>The instruments may not be used for any other than the described indication or application area.</p> <p>Excessive temperatures due to insufficient water cooling must be avoided (possible damage of pulp)</p> <p>The indicated speed may not be exceeded (risk of fracture/injury)</p> <p>Jamming or using the instrument as a lever must be avoided (risk of fracture/injury)</p> <p>Processing of soft materials must be avoided</p>
<b>Shelf Life</b>	10 years
<b>Number of Reworks</b>	The recommended number of reworks is a maximum of 15 times.
<b>Sterilization</b>	Products are not Sterile

## 2. WARNINGS

- Be sure to read the user manual before using.
- Avoid any contact with H<sub>2</sub>O<sub>2</sub> (hydrogen peroxide).
- Make sure that specified immersion times in the cleaning and disinfecting agents are not exceeded.
- Operative burs must be sterilized before each use to avoid potential transfer of infectious organisms.
- Check whether there is any moisture on the packaging after sterilization. Use the product in packaging that does not cause moisture.
- Store products out of children's reach.
- For dental use only.
- Dispose of as contaminated medical waste.
- Any serious adverse event occurring during use and/or life cycle of the device should be reported to the manufacturer, notified body and the local component authority.

## 3. METHOD OF USAGE

### 3.1 Proper Use

- The instruments have to be prepared prior to first use.
- Make sure that only technically and hygienically perfect and cleaned power systems (hand pieces and contra-angles) are used.
- Depending on their shank type, insert the instrument into the chuck as deeply as possible. Make sure that they are properly locked.
- The instrument must be rotating at the desired speed before contact is made with the tissue or material.
- Avoid jamming and using the instrument as a lever. Excessive contact pressure has to be avoided.
- For eye protection wear safety glasses. Use appropriate respiratory protection (mouth and nose). In the dental laboratory, use appropriate suction unit.
- Avoid contact with the instruments' working parts as this may increase the risk of injury.

Make sure to provide sufficient cooling by means of air/water spray. Additional external cooling is required when using instruments with extra-long shanks or oversized working parts. Damaged or corroded instruments have to be discarded.

### 3.2 Recommended speeds

HEAD DIAMETER (ISO $\varnothing$ (1/10 mm))	FG(Max. RPM)	HP(Max. RPM)
007-014	300,000	160,000
016-023	280,000	120,000
025-045	180,000	100,000
047-065	80,000	80,000
066-093	60,000	50,000

For optimum results when using Bosphorus branded burs, 50% of the maximum speed should be used.

### **3.3 Contact Pressure**

Excessive contact pressure (> 2N) has to be avoided.

In cutting instruments, this can lead to damage to the working part and to chips of the diamond as well as an excessive generation of heat.

### **3.4 Reprocessing**

Those instruments that are supplied non-sterile have to be prepared prior to first use.

The service life of tools is not unlimited. The service life of instruments is normally determined by wear and tear from their intended surgical use, not by their rework. Instruments should be meticulously inspected prior to each use to ensure they are fully functional.

## **4. PRODUCT CLEANING**

### **4.1. Required materials and equipment**

Washers-Disinfectors (WDs) according to ISO 15883-Cleaning Detergents

Detergents with a pH between 7-10 are recommended; The maximum pH should be 12.

### **4.2. Manual cleaning and disinfections**

Make sure that the burs do not come into contact with each other during this process.

#### **Pre-Cleaning**

1. Rinse the contaminated device under cold running tap water for at least three minutes. Use a soft brush to help remove most dirt and debris.

#### **Manuel cleaning**

2. Immerse the device in an alkali cleaner or detergent solution for at least 10 minutes. Follow the enzymatic cleaner or detergent manufacturer's instructions for use for the correct exposure time, temperature, water quality, and degree of concentration.
3. Manually clean the device until no residue is visible on its surface. Use a soft brush to remove dirt and debris. If applicable, activate joints and all other movable device features to expose all areas to the detergent solution.
4. Prepare an ultrasonic bath with a fresh detergent solution using an enzymatic cleaner or detergent at the concentration and temperature specified in the detergent manufacturer's instructions.
5. Ultrasonically clean instruments using a minimum frequency of 30 kHz for at least 10 minutes.
6. Visually check the device. Repeat steps 2–5 until no visible residue remains on the device.

#### **Manuel disinfection**

7. Immerse the device in an enzymatic cleaner/ disinfection solution for at least 15 minutes. Follow the manufacturer's directions for use for the correct exposure time, temperature, water quality, and degree of concentration.
8. Rinse the device thoroughly in purified water for at least two minutes.
9. Dry the instruments.

### **4.3. Automatic cleaning / disinfection**

#### **Minimum program parameters**

1. Pre-rinse with cold tap water for at least 3 minutes
2. Main wash at least 2 minutes below 40°C. add detergent
3. Washing with detergent at 55°C for at least 5 minutes
4. Pre-rinse for at least 1 minute with tap water above 40°C
5. Rinsing with cold demineralized water for at least 1 minute
6. Thermal disinfection with demineralized water above 90°C for at least 5 minutes
7. Drying above 90°C for at least 40 minutes
8. During evacuation, check that the tools are completely free of visible dirt. If necessary, repeat the program or perform manual cleaning.

### **4.4. Drying**

Dry the burs for 5 minutes with paper towels or using dry heat not exceeding 140°C.

### **4.5. Controls**

After each cleaning/disinfection, the products must be completely free of visible residues and other contaminants. After cleaning, inspect the burs, preferably with a magnifying glass. If necessary, the cleaning process should be repeated until the product is visually clean.

**Function:** Before each use, it should be checked for cracks, cracks, deformation, corrosion and proper function.

Check for possible damage; Do not use drills that are oxidized, worn, or whose diamond coating has peeled off.

### **4.6. Packaging and Labeling**

Package products in sterile barrier system packaging suitable for steam sterilization or in a sterilization packaging or equivalent container as soon as possible after reprocessing in accordance with ISO 11607 and label the packaging in accordance with ISO 15223-1.

### **4.7. Sterilization**

The sterilization process must be carried out in accordance with ISO 17665 and steam sterilizers operating under vacuum and kept under the required conditions.

We recommend using the following steam sterilization parameters:

- **Minimum holding time of 5 minutes at 134°C**

## 5. POTENTIAL EFFECTS OF USING WORN INSTRUMENTS

- Discard any instruments with damaged or worn working parts as the use of damaged or worn instruments would have a negative effect on the work result.
- Damaged and deformed burs will cause the instrument to vibrate and lead to poor preparation margins and rough surfaces.
- Please avoid excessive contact pressure during use as this may result in damage to the working part (Nicks on the blades, premature blunting and excessive heat generation).
- To avoid undesirable heat generation during preparation, make sure to provide sufficient cooling by means of air/water spray (at least 50 ml/min).
- Additional external cooling is required when using instruments with a total length of more than 22 mm or a head diameter exceeding 2,5 mm.
- Not observing the maximum permissible speed will result in an increased risk of injury.
- There is an increased risk of infection in reusable instruments which have not been properly reprocessed.

## 6. STORAGE

- Tools should be stored in a clean, dry environment protected from direct sunlight, germs/insects, extreme temperatures and humidity.
- Do not store in the same room as solvents or chemicals.
- Do not damage the packaging material.
- The product is subject to improvement/change without notice. For this reason, it is recommended to use the first purchased products first.

## 7. SAFETY AND POSSIBLE SIDE EFFECTS

The above mentioned recommendations with respect to cooling, contact pressure, disinfection, cleaning and sterilization are to be strictly observed. The instruments should only be used for the intended application, as per the symbolic identification.

Non-observance of these safety recommendations may lead to damage of the power system and/or injury, such as thermal necrosis, undesired preparation of tissue, damage to tissue or nerves, violation of the biological width, or infections.

### MRI Safety



Diamond dental burs are not implantable devices.

## 8. LIABILITY

It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Bosphorus partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

## 9. SYMBOLS OF LABEL

	5.7.7 Medical device		5.4.4 CAUTION
	5.7.10 Unique device identifier		5.2.7 Non-Sterile
	5.1.1 Manufacturer		5.3.1 Fragile, handle with care
	5.1.6 Catalogue Number		5.3.4 Keep dry
	5.1.5 BATCH CODE		5.3.2 Keep away from sunlight
	5.1.3 Date of manufacture		5.2.8 Do not use if package is damaged and consult instructions for use
	5.1.11 Country of manufacture	<b>ISO NO</b>	According to EN ISO 6360-1,2,3,5,6,7
	5.4.3 Consult Instructions For Use Electronic Instructions For Use	<b>TYPE</b>	Defines type of operative burs
	5.1.4 Expiration date		